



High flow pump, two split system



High flow pump, two split system



Pros:

Flow split close to the trap ensures fast gradient transitions as the flow rate is 200-300 $\mu\text{l}/\text{min}$ up to the first split.




Cons:

The flow rate changes from loading to analysis may cause pressure fluctuations affecting the packing of the trap/column or induce leaks

If the first split leaks during loading you'll lose sample

1. Optional in-line solvent filter (Upchurch A314 with 2 μm peek frit A702) connected via peek tubing (127 μm ID) to reduce risk of clogging downstream lines/columns
2. Sample loop; e.g. PEEKsil tubing 15 cm x 1/16" x 0.3 mm ID: 10.603 μl (Upchurch part#630015)
3. Injection needle (home made for Spark Holland Endurance AS 100 μm ID x 37 cm: 3 μl)
4. Transfer line fused silica 50-100 μm ID x 25 cm: 0.5-2 μl
5. Peek MicroTee (Upchurch P-775 or P-875 w/ mounting whole)
6. Trap column: e.g. fused silica 100 μm ID x 20 cm = 1.6 μl (PicoTip Integrafrit # IF360-100-50-N-5) packed with MagicC18AQ 200 \AA 5 μm c.a. 2-4 cm long
7. Peek MicroCross (Upchurch P-777), high voltage applied through 0.5 platinum or gold wire
8. Empty tip or separation column: e.g. fused silica 75 μm ID x 10-60 cm tip pulled manually with microflame torch, packing MagicC18AQ 100A 5 μ 10 cm long
9. Flow split : fused silica 25-50 μm ID x 15-30 cm open in detector position; adjust ID and length to regulate flow rate through column to 200-300 nl/min
10. Flow split : fused silica 100 μm ID x 15 cm open in waste position

High flow pump, single constant open split system



High flow pump, single constant open split system



Pros:

- Constant flow rate at pump (200-300 $\mu\text{l}/\text{min}$)
- Reduced the risk of pressure fluctuations
- Self regulated flow rate through trap and column
- Reduced risk of sample loss during loading




Cons:

- Increased void volume leads to increased delay time

1. Optional in-line solvent filter (Upchurch A314 with 2 μm peek frit A702) connected via peek tubing (127 μm ID) to reduce risk of clogging downstream lines/columns
2. Peek MicroTee (Upchurch P-890), NOTE: mount as close to the AS valve as possible to minimize void volume, use small ID line to connect to AS valve (e.g. 5cm x 127 μm ID = 630 nl, 5cm x 50 μm ID = 98 nl)
3. Flow split : PEEK or fused silica 25-50 μm ID x 15-30 cm open in detector position; adjust ID and length to regulate flow rate through column to 200-300 nl/min
4. Sample loop; e.g. PEEKsil tubing 15 cm x 1/16" x 0.3 mm ID: 10.603 μl (Upchurch part#630015)
5. Injection needle (home made for Spark Holland Endurance AS 100 μm ID x 37 cm: 3 μl)
6. Transfer line fused silica 50-100 μm ID x 25 cm: 0.5-2 μl
7. Peek MicroTee (Upchurch P-775 or P-875 w/ mounting whole)
8. Trap column: e.g. fused silica 100 μm ID x 20 cm = 1.6 μl (PicoTip Integrafrit # IF360-100-50-N-5) packed with MagicC18AQ 200 \AA 5 μm c.a. 2-4 cm long
9. Peek MicroCross (Upchurch P-777), high voltage applied through 0.5 platinum or gold wire
10. Empty tip or separation column: e.g. fused silica 75 μm ID x 10-60 cm tip pulled manually with microflame torch, packing MagicC18AQ 100A 5 μ 10 cm long
11. Flow split : fused silica 100 μm ID x 15 cm open in waste position



e.g. nanoAcquity





Pros:

- Low flow, less waste!
- Improved peak capacity and peak shape
- High chromatographic reproducibility



Cons:

- Flow rate needs to be adjusted for every new column
- More sensitive to solvent impurities
- more difficult to find leaks

1. Sample loop (ss): 5 or 10 μ l
2. Transfer line fused silica 25-40 μ m ID x 25 cm: 0.1-0.3 μ l
3. Peek MicroTee (Upchurch P-775) closed with one plug
4. Trap column: e.g. fused silica 100 μ m ID x 15 cm = 1.18 μ l (PicoTip Integrafrit # IF360-100-50-N-5) packed with MagicC18AQ 200 \AA 5 μ c.a. 2-4 cm long (NOTE we reuse the Integrafrit by flushing the beads out using the HPLC)
5. Peek MicroCross (Upchurch P-777), high voltage applied through platinum or gold wire
6. Empty tip or separation column: e.g. fused silica 75 μ m ID x 10-60 cm, tip pulled manually with microflame or laser puller, packed with MagicC18AQ 100 \AA 5 μ 10-60 cm long (avoid any void volume between trap and column, cut the column to desired length, such that beads are packed all the way to the end of the fused silica, most commonly used column length at UWPR is 20-30cm)