NSI Source: 3-D printed carriers

Advancing Proteomics

4-40


6-32

6-32

4-40

UWPR

When using stainless steel crosses of tees there is the potential for the high voltage to arc. To minimize that risk we use 3D printed carriers, using PETG filament to print the carrier.



Carrier with P-777 cross connected to high voltage lead and 360µm OD column, Trap column and waste line

Optional, to hold PEEK tee with mounting hole, use 4-40 screw 18-8 Stainless Steel Female Hex Thread Adapter, 4-40 to 6-32 Thread Size (98434A102, McMaster)

> Nylon Raised-Head Thumb Screws 8-32 Thread Size, 1" Long, black (94323A598, McMaster)

Small Adjustable Clamping Arm, 6-32 threaded Post (PM3, Thorlabs)


18-8 Stainless Steel Male-Female Hex Thread Adapter, 6-32 to 4-40 Thread Size (92499A037 McMaster)

> 3D printed carrier for Idex PEEK cross (P-777) or PEEK tee (P-875) Or VHP cross (UH-752) or tee (UH-750 or UH-753)



98434A102 18-8 Stainless Steel Female Hex Thread Adapter. 4-40 to 6-32 Thread Size (McMasterCarr) 94323A598 Nylon Raised-Head Thumb Screws 8-32 Thread Size, 1" Long, black PM3 Small Adjustable Clamping Arm, 6-32 threaded Post (Thorlabs) (McMasterCarr) 92499A037 18-8 Stainless Steel Male-Female Hex Thread Adapter, 6-32 to 4-40 Thread Size 3D printed carrier for Idex PEEK cross (P777) or PEEK tee (see separate files to print this carrier) Rail Carrier, Counterbored Hole 1"x 1" (Thorlabs)

(McMasterCarr)



from below to hold hex thread adapter to the RC1 rail carrier

University of Washington Proteomics Resource